BLANCONE[®]

ARCUS⁺

MOBILE

LED Bleaching Lamp


Lampada sbiancante a LED

LED Bleaching Lampe

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LED Bleaching lamp MOBILE

GB

I.	Unit description and functions	2
II.	Technical specifications.....	2
III.	Packaging / components	3
IV.	Safety precautions	3
V.	Installation.....	6
VI.	Preparation for operation	8
VII.	Daily care and maintenance	9
VIII.	Preparation for work with patients.....	10
IX.	Positions of bleaching head	11
	against patient's teeth	
X.	Determination of distance	12
	to patient's teeth	
XI.	Problems and solutions	13
XII.	Warranty	14
XIII.	Symbols	85
XIV.	Service data	85
XV.	Declaration of conformity	86
XVI.	Led bleaching unit data	87

Caution!

Before installing and operating with the unit, carefully read this manual.

I. UNIT DESCRIPTION AND FUNCTIONS

BLANCONE ARCUS bleaching device is designed to bleach the solid tooth tissue by photoactivation of a whitening gel applied in dentistry. The appliance should only be used by a dental practitioner and in a dental practice. The device consists of a Bleaching Head, Stand and a Power Adapter. BLANCONE ARCUS BUILT-IN is manufactured in accordance with the requirements of Directive 93/42/EEC.


II. TECHNICAL SPECIFICATIONS

1. Operating voltage to the power adapter - 100-240V / 50-60Hz, to the whitening head - 24VDC.
2. Current consumption - 0.4-0.2 A
3. Dimensions:
 - Bleaching head
 - length - 220 mm
 - width - 155 mm
 - height - 95 mm
 - Stand
 - length (assembled) - 1300 mm
4. Weight:
 - Bleaching head - 640 g
 - Stand - 8500 g
5. Irradiance - over 250 mW / sq.cm. measured in the output window.
6. Operation time - 10 sec. - 30 min / $\pm 5\%$ /
7. Emitted light - blue 430 - 490 nm.
8. Light source - 10 LED modules x 5W = 50W
9. Fan air cooling.
10. Ability to start and stop the Bleaching Head at any time.
11. Maximum continuous operation time 99 min – after every 99 minutes of continuous operation let the unit cool for at least 10 minutes.

Upon request the manufacturer of this device will provide all other necessary documentation and/or information for the user's technical staff to repair those parts of the apparatus which the manufacturer has indicated as suitable for repair outside the manufacturer's aftermarket service facility.

III. PACKAGING / COMPONENTS

1. Power cord
2. Power Adapter 100-240VAC / 24VDC
3. Horizontal arm
4. Central position locking mechanism
5. Front position locking mechanism
6. Base with 5 wheels
7. Vertical arm
8. Bleaching head
9. Counterweight
10. Control panel
11. "+" Timer or "Pause" button
12. "-" Timer or "Pause" button
13. Start / Stop button
14. Display
15. Bleaching head
16. Cables and connectors
17. Safety glasses - 2 pcs.

The device should only be moved when all 5 wheels are attached to the Base 6.

Positioning of the horizontal arm is done by the Central position locking mechanism 4.

Fine positioning is done by turning the Bleaching head so that the light is perpendicular to the patient's teeth and centered on them. If the movement of the Bleaching head is difficult or it does not hold the selected position it is necessary to adjust the Front position locking mechanism 5.

IV. SAFETY PRECAUTIONS

1. The device must be used in accordance with the User Manual. It is not allowed to use it for any other purpose or by unauthorized personnel.
2. The device is equipped with European type plug for voltage 100-240V/50-60Hz. Usage of the device with other plugs and/or other voltages can damage it or may lead to certain malfunctions.

Electrical safety is provided by Class I protection against electric shock according to EN 60601-1.

3. BLANCONE ARCUS should be operated indoors under the following conditions:
 - temperature from + 10 ° to + 40 ° C;
 - relative humidity of the air - 30 - 75%; The device should not be used in rooms with high humidity levels, condensation or wet floor, walls and/or ceiling.
 - atmospheric pressure 700 - 1060 hPa;
 - absence of chemically active and flammable substances.
4. Strong electromagnetic fields in the building may cause interference and malfunction. If their source cannot be determined, change the location and turn it on in another electrical plug.
5. BLANCONE ARCUS is a source of extremely intense light in the blue light spectrum to which the human eye has a high sensitivity. Serious measures need to be taken in order to protect patients, medical staff and accidentally close people, animals and plants. Always use protective glasses for the operator and for the patient's eyes, wear mask and facial sunscreen. Eye irritation carries a high risk. Never direct the light to the eyes! Irradiation must be limited to the place of work.

Caution!

The light emitted is at a wavelength of 430-490 nanometers.

Use safety glasses that comply with the following requirements:

- tightly cover the eyes and the temple even if the operator bears dioptric glasses.
 - made from a voluminous colored impact-resistant plastic.
 - admit light with a wavelength exceeding 600nm.
 - reduce the intensity of the blue spectrum over 100 times.
 - have a stable mechanical structure and do not have any scratches or damage on their surface.
6. The device should not be used with pregnant and breast-feeding mothers, patients with: oncological problems, severe parodontal pathology, recessions, dental paresthesia and under 17 years of age.
The device may only be used after medical consultation on or by persons with implanted cardiac pacemaker; persons suffering from photobiological reactions; taking photosensitive medications; individuals with cataract surgery, people with retinal illnesses, allergy sufferers; persons with multiple sensitive skin or dermatitis, etc.
If the bleaching protocol is not followed, pain, hypersensitivity, enamel defects and even burning of the non-calcified, soft tissue may occur.
 7. It is desirable to check the sensitivity of the patient to the intensity of the light. The device has two levels of intensity (High and Low) and the optimal conditions are selected by adjusting the distance of the device to the patient's teeth. Testing should begin at a high level "Hi" (High) and a distance of 5 cm. If the thermal sensation is unpleasant and irritating, increase the distance to achieve an acceptable thermal comfort. If the thermal sensation is still unpleasant at 10-15 cm switch to low intensity - "Lo" (Low) and shorten the

distance to the teeth. If the patient has no sensitivity or sensitivity cannot be measured for some reason, longer distance should be used.

8. Do not be placed any solvents, flammable liquids and heat sources near the device and/ or the power adapter, because exposing it to any of these may cause damage, fire or an explosion.
9. Cleaning agents should not penetrate the device, as this may lead to malfunctions, damage or fire.
10. Protect the outer insulation of all cables from damage that may be caused by sharp objects, strong pulling, rodents and chemical substances. If any such damages are noticed switch off the device and immediately take it to the service facility.
11. In the event of a thunderstorm, the device must be stopped and disconnected from the power network.
12. Do not force the horizontal arm to avoid mechanical shocks to the counterweight or the bleaching head. In the event of a mechanical shock, if the bleaching head is damaged, the device should not be used. It must be taken immediately to the service facility.
13. The movement of the bleaching head and horizontal arm without adjustments to the Central and/ or Front position locking mechanism may damage the device.
14. Any repairs may only be performed by service technicians authorized by the manufacturer.
15. When replacing damaged components, only original BLANCONE ARCUS spare parts should be used. The warranty of the device does not cover damages caused by the use of non-original spare parts. The device or any of its components should not be disassembled while plugged into the power network!
16. Before each procedure, the bleaching head window should be disinfected with disinfecting solution.
17. Take precautions against accidental injury to the human body from mechanical parts (moving and stationary), pinching, inertial acceleration of the counterweight etc. Safe handling of the device is vital to avoid injury.
18. Before the procedure each patient must sign an informed consent, detailing the possible side effects of the process. Note: It is also advisable to visit a suitable Teeth whitening in dental practice course.
19.  Fragile! Use caution when transporting, using and storing the device! When transported in an assembled condition, two persons are required to keep the moving parts. With any unevenness on the flooring, the device should only be moved by two people.
20.  According to Directive 2012/19/EEC, this symbol indicates that the product should not be disposed as a general waste at the end of its lifespan. The product must be taken to a specialized center for the separate collection of electrical and electronic equipment. Proper disposal of equipment that is no longer used prevents negative consequences for the environment and human health!

V. INSTALLATION

BLANCONE ARCUS is intended for use as an autonomous appliance. Use the device indoors and only on a horizontal surface (floor).

1. Unpack the transport boxes, remove the components and initially assemble the base with the 5 wheels.

2. Mount the vertical arm to the center of the base.

3. Mount the horizontal arm of the stand and fix it with the side stop screw using a 3 mm Allen key (included in the package).

4. Place the counterweight at the rear of the horizontal arm by rotating it clockwise. Note: It is important to mount the counterweight securely on the arm (until it stops rotating clockwise).

5. Mount the bleaching head at the front of the horizontal arm in the following order:
 - Loosen the adjustment ball (rotate counterclockwise)
 - Tighten the bleaching head until it stops.
 - Tighten the adjustment ball (rotate clockwise) so that the bleaching head can move in all directions and stay in each set position.
 - Tighten hard the adjustment ball to fix the bleaching head.


- 6.** The horizontal arm allows movement in vertical direction. A fixing ball locks the horizontal arm and the bleaching head into a suitable position.

- 7.** Connect the cable from the bleaching head to the one from the front of the horizontal arm using the connectors (see figure). Note: the device will not work if the cable connectors are incorrectly coupled.

- 8.** Connect the cable from the back of the horizontal arm to the one on top of the vertical arm using the connectors (see figure). Note: the device will not work if the cable connectors are incorrectly coupled.

- 9.** Connect the cable on the bottom of the vertical arm with the power adapter using the connectors (see figure). Note: the device will not work if the cable connectors are incorrectly coupled.

- 10.** Plug the power cord into the power network.

VI. PREPARATION FOR OPERATION

1. The BLANCONE ARCUS bleaching device should be placed in a suitable position around the patient and the dental unit. See the positioning examples (p.IX). The bleaching head must be in the same vertical plane and parallel with the arc of the teeth with a distance between 5-15 cm. See the positioning examples (p.X).
2. Plug the power adapter into the power network. Only one point/dot on the display lights up. Press START / STOP button and then the "Hi" (High) or "Lo" (Low) intensity indication appears.
3. The display shows "Hi" - a high intensity level. If you want to switch to the "Lo" (Low) intensity mode, press the TIMER (+) or TIMER (-) button once. Then the display will show "Lo". Every time before setting the required time, you will be able to select the Hi or Lo intensity level.
4. Press the START / STOP button once again and set time option will be displayed. Note that the device will save the time frame used in the latest bleaching procedure. Required time frame for operation is set by pressing TIMER + (increasing time) and TIMER- (reducing time):
The display on the bleaching head will indicate timeframes between 10 sec. and 9 min. and 50 sec. with a comma between the two digits.
Examples:
Display indicates: 13. This means the timeframe is set to 13 min.
Display indicates: 0.1. This means the timeframe is set to 10 sec.
Display indicates: 9.5. This means the timeframe is set to 9 min. and 50 sec.
The display on the bleaching head will indicate timeframes between 10 min. and 30 min. without a comma between the two digits.
5. Press the START / STOP button one final time and the device will start to emit light. A flashing dot on the display will indicate that the device is working properly.
6. If either of the TIMER + / TIMER - buttons is pressed during operation, the unit switches to Pause mode. A blinking display with the remaining time will indicate that the device is paused. Operation is resumed by pressing the TIMER + / TIMER button again.

7. To Stop the bleaching device at any time press the START / STOP button. If the device is not manually stopped it will stop after the selected time has elapsed. Note that the fans will continue to work for some time in order to cool down the unit.
8. At the end of the day turn off the device by holding down the Start / Stop button for 3-4 seconds or by disconnecting it from the power network.
9. In case of overheating the thermal protection is triggered and the device will cease to emit light.
The display will indicate "Oh" and the fans will turn on for 1 minute.
After cooling, the device will continue its normal operation.
Overheating can only happen in an emergency situation or by accidental clogging of the cooling vents.
That's why the cooling vents located on the bottom of the bleaching head must always be unobstructed in order to allow air circulation.

VII. DAILY CARE AND MAINTENANCE

1. Disinfect the PVC window on the bleaching head before and after each patient with cotton or other nonabrasive cloth soaked in alcohol.
2. To disinfect the device and its parts, spray disinfectant onto a piece of soft cloth or cotton and clean it.
Do not use abrasives or solvents as this may cause damage.


VIII. PREPARATION FOR WORK WITH PATIENTS

1. Before starting work with the bleaching device, insulate the patient's non-calcified soft tissue, place a cover or protective napkin on the face, apply protective UV cream and put on protective glasses. The bleaching gel is applied according to the directions for use outlined by the gel manufacturer. The patient should be under constant supervision by the dentist during the bleaching process in order to avoid discomfort and to ensure all safety and technological procedures are applied.
2. It is desirable to check the sensitivity of the patient to the intensity of the light. The device has two levels of intensity (High and Low) and the optimal conditions are selected by adjusting the distance of the device to the patient's teeth. Testing should begin at a high level "Hi" (High) and a distance of 5 cm. If the thermal sensation is unpleasant and irritating, increase the distance to achieve an acceptable thermal comfort. If the thermal sensation is still unpleasant at 10-15 cm switch to low intensity- "Lo" (Low) and shorten the distance to the teeth. If the patient has no sensitivity or sensitivity cannot be measured for some reason, longer distance should be used.



IX. POSITIONS OF BLEACHING HEAD AGAINST PATIENT'S TEETH

CORRECT


INCORRECT


X. DETERMINATION OF DISTANCE TO PATIENT'S TEETH

HIGH intensity

5 cm distance


9 cm distance


15 cm distance

LOW intensity

5 cm distance

9 cm distance

15 cm distance

Hot

Warm

Acceptable

Weak

XI. PROBLEMS AND SOLUTIONS

The device does not work.

Check if the power cord from the adapter is plugged into the power network and all cable connectors are correctly coupled.

Damaged power cord.

Do not use the device! Switch it off and immediately take it to the service facility

Unpleasant or irritating thermal sensation by the patient.

Increase the distance of the device to the patient's teeth to 10-15 cm. If the thermal sensation is still unpleasant at this distance switch to low intensity- "Lo" (Low) and shorten the distance to the teeth.

The display is blinking and no light is emitted.

The device is in PAUSE mode. To resume operation, press either + or - buttons on Bleaching Head.

The Bleaching Head cannot be fixed.

Read Point 5. In Section V. INSTALLATION and slightly adjust the tension of the Bleaching Head until the desirable mobility is reached.

The display only shows "Oh".

In case of overheating the thermal protection is triggered and the device will cease to emit light. The display will indicate "Oh" and the fans will turn on for 1 minute. After cooling, the device will continue its normal operation. Overheating can only happen in an emergency situation or by accidental clogging of the cooling vents. That's why the cooling vents located on the bottom of the bleaching head must always be unobstructed in order to allow air circulation.

For all other questions regarding the installation and/or operation of the device, please contact the manufacturer or your local authorized representative.

XII. WARRANTY

1. The warranty period of the BLANCONE ARCUS BUILT-IN bleaching device is 24 (twenty-four) months from the date of purchase. If the purchase date is not filled in Section X. LED BLEACHUNG UNIT DATA, the warranty begins on the date of production.
2. During the warranty period, the replacement of the defective elements shall be done free of charge by the manufacturer. Note - LED modules have a 6-month warranty.
3. The device and its components must only be used as intended by the manufacturer purpose outlined in this User Manual. Any other use of the device will lead to a warranty lose and no obligations from the manufacturer.
4. If the device fails as a result of improper operation during the warranty period (mechanical, chemical, thermal, electrical), non-intended use, inappropriate storage etc., warranty will be void and the repair cost will be at the expense of the user.

The device should not be used with damaged cables. If any such damages are noticed switch off the device and immediately take it to the service facility.

If liquids such as water and solvents, aggressive or flammable substances, vapors, insect or rodent enter the device, switch it off and immediately take it to the service facility. Damages from the above will void the warranty.

No claims shall be accepted for damages arising from electrical shocks, thunderstorms, non-compliance with electrotechnical safety measures or insufficient protection of patients, staff, other people, animals, plants, and objects from light radiation.

The manufacturer does not owe any compensation for lost profits during the period when the device is damaged or malfunctioned, no matter what the reason.

Claims for damages and claims arising out of non - compliance with the bleaching procedure are not accepted. Including but not limited to: shorter or longer time period use of the device with the bleaching material set by the gel manufacturer; unsatisfactory bleaching result; patient's harms from the bleaching material; overdosage; insufficient insulation of patient's non-calcified soft tissue; improper protection of the patient and staff; inappropriate, expired, intended for another wavelength and/or with unsuitable concentration bleaching material.

The warranty is forfeited and no claims shall be accepted for damages as a result of incorrect or insufficient care and attention for protection during transport, unpacking, moving, handling and storage of the device.

In the event of disputes arising out of the application and interpretation of this User Manual, those will be settled by the courts in the city of Plovdiv, under the current Bulgarian legislation.

5. The warranty of the device will be forfeited if any repairs or modifications are done by unauthorized personnel outside of the manufacturer's service facilities and/or non-original spare parts are used.
6. The manufacturer recommends that customers should check whether the device parameters are within the permissible limits once a year. Technical condition tests and verifications may only be carried out at the manufacturer's service facilities or by an authorized representative.
7. All repairs must be done at the manufacturer's service facilities with the following address:

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XIII. SYMBOLS

SIMBOLI - SYMBOLE - SYMBOLES - SÍMBOLOS - SÍMBOLOS

	Manufacturer BG LIGHT LTD, 155, Vasil Aprilov blvd., 4027 Plovdiv, Bulgaria
	Pursuant to Directive 2012/19/EU, this symbol shows that the product must not be disposed of as urban waste at the end of its operating life.
	Pay attention to the instructions accompanied by this symbol.
	Fragile!
	Product with CE mark
	Applied part type B
Ω	Ohm (electric resistance unit)
s	Second (time unit)
W	Watt (power unit)
Hz	Hertz (frequency unit)
mm	Millimeter (length unit)
A	Amper (electric current unit)
V	Volt (electric voltage unit)
Pa	Pascal (atmosphere pressure unit)
°C	Degrees Centigrade (temperature unit)
g	Gram (weight unit)
SN	Serial number
LOT	Lot

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Product: **BlancOne ARCUS* Mobile**
Classification: Active device in Class IIa of the Medical Device Directive 93/42/EEC+AC:2007

Classification is done by the manufacturer according to Medical Device Directive 93/42/EEC+AC:2007, Rule 9, Annex IX.

Notified body: TUV NORD Polska Sp. z o.o., ul.Mickiewicza 29, 40-085 Katowice, Poland.

This Declaration of conformity is valid only in combination with our certificates of Notified body TUV NORD Polska Sp. z o.o.
Certificates N: AC090 100/0951/4047/2014, AC090 MD/0951/4047/2014, TNP/MDD/0122/4047/2015..

The manufacturer declares under its sole responsibility that the products are developed and produced in conformity with MDD 93/42/EEC+AC:2007 and the following applicable standards:

EN 60601-1	Electrical safety,
EN 60601-1-2	Electromagnetic compatibility
EN 62304	Medical device software. Software life cycle processes.
EN 62353	Medical electrical equipment. Recurrent test and test after repair of medical electrical equipment
EN 62366	Medical devices. Application of usability engineering to medical devices.
EN ISO 14155-1	Clinical investigation of medical devices for human subjects. General requirements.
EN ISO 14155-2	Clinical investigation of medical devices for human subjects. Clinical investigation plans.
EN ISO 14971	Medical devices – Application of risk management to medical devices
EN ISO 15223-1	Medical devices -- Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
EN 1041	Information supplied by the manufacturer with medical devices
EN ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process

All company products are manufactured under the current Quality Management System, ISO 9001 and ISO 13485.

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Date of production:

Date of purchase:

**The quick and gentle
professional teeth whitening**

www.blancone.eu/pro

blancone

BlancOne Team PRO

blancone_ids

blancone

IDS Spa

17100 Savona - Via Valletta San Cristoforo, 28/10

Tel. +39 019 862080 - Fax +39 019 2304865

www.ids dental.it - info@idsdental.it

**The quick and gentle
professional teeth whitening**

www.blancone.eu/pro

IDS Spa

17100 Savona - Via Valletta San Cristoforo, 28/10
Tel. +39 019 862080 - Fax +39 019 2304865
www.idsdental.it - info@idsdental.it